



Less Mass, More Science

How cavity-enhanced spectroscopy can improve science returns and decrease cost for in-situ planetary measurements

Bradley M Gibson, PhD, 389R

© 2017 California Institute of Technology. Government sponsorship acknowledged.




Jet Propulsion Laboratory
California Institute of Technology

- What is Tunable Laser Spectroscopy?
 - What can we measure with TLS?
 - What goes into making a TLS for planetary missions?
 - How can we improve, and what could we do then?
-
- What is cavity-enhanced spectroscopy?
 - What are the design tradeoffs?
 - What is unique about our design?
 - How well does the spectrometer perform?
 - What can we improve in the future?

Tunable Laser Spectroscopy Overview

What is TLS?

Tunable Laser Spectroscopy (TLS) uses the attenuation of narrow-band laser light to achieve **sensitive, selective** detection of molecules over a **wide range of pressures**.



What can we measure?

We can measure molecules that are:

- Between 2 and 12 atoms in size (more if rigid, like benzene)
- Present at parts-per-billion volume or greater (ppb wt for solids)
- In the gas phase, or can be vaporized

TLS can:

- Easily distinguish species of similar mass – e.g. $\text{C}^{13}\text{O}^{16}\text{O}^{16}$ from $\text{C}^{12}\text{O}^{16}\text{O}^{17}$ and C^{12}H_4 from O^{16}
- Provide highly localized measurements

What are the real-world applications?

Earth Science:

- Observe water transport through tropopause / stratosphere via isotopic measurements
- Measure methane / CO₂ isotope fluxes to localize and quantify sources and sinks
- Monitor concentrations of reactive species (halogen compounds, NO, OH) and ozone
- Monitor emission of volcanic gases

What are the real-world applications?

Planetary Science:


- Isotope ratio studies of various compounds – CO_2 , H_2O , CH_4 , PH_3 , NH_3 , SO_2 – as tracers of solar system formation, geological and atmospheric processes of each body
- Localized mixing ratios for atmospheric layers or near sites of geological or atmospheric activity
- Detection of trace quantities of biologically-relevant gases in atmosphere or processed samples

What are the real-world applications?

Miscellaneous:

- Laboratory measurements in support of astronomy
- Observations of combustion chemistry / physics
- Breath analysis for medical diagnostics
- Safety monitoring for mining, drilling, and manufacturing industries

What is TLS?



$$T = \frac{I}{I_0} = 50\%$$

$$A = \log_{10} \frac{I_0}{I} = 0.3$$

What is TLS?

$$A = \sigma \rho l$$

Absorption cross-section (cm^2)

Path length (cm)

Volume density (cm^{-3})

$$\alpha \equiv \sigma \rho$$


Absorption coefficient (cm^{-1})

Noise-equivalent absorption (NEA):



Absorption signal equal to one standard deviation of the absorption measurement.

- Detector and shot noise
- Laser wavelength and intensity fluctuations
- Mechanical and thermal perturbations
- Etalons / fringing



What is TLS?




What is TLS?




What is TLS?




What is TLS?




What is TLS?




What is TLS?




What does a TLS system look like?




What does a TLS system look like?





What does a TLS system look like?




What does a TLS system look like?





Example: SAM suite on Curiosity



Example: SAM suite on Curiosity



Example: SAM suite on Curiosity



What can we improve?

We can increase path length.

We can decrease noise.

We can decrease size / weight of instrument.

We can decrease sample volume / gas handling requirements.

What new science goals can we enable?

We can measure low-abundance targets without pre-concentration, such as methane on Mars or clumped isotopes.

We can measure difficult-to-collect samples, such as water vapor from comets or plumes on Enceladus.

We can operate in atmospheres containing difficult-to-pump balance gases, such as helium on Saturn.

Cavity Enhanced Spectroscopy

What is cavity-enhanced spectroscopy?

Cavity-enhanced spectroscopy uses an optical resonator to let light traverse the same path thousands of times, dramatically increasing path length.




What is cavity-enhanced spectroscopy?




$$\text{EPL} = \frac{L}{1 - R}$$

$$\frac{12 \text{ cm}}{1 - 0.99998} = 600 \text{ m}$$




Cavity resonance conditions




Resonance condition:
$$L = \frac{n\lambda}{2}$$


Cavity resonance conditions



What is cavity-enhanced spectroscopy?



What is cavity-enhanced spectroscopy?




Off-Axis Injection:



Injecting light at an offset from the center of the mirrors increases distance traveled before beam overlaps (reentrance).

- Decreases FSR, noise from inconsistent coupling
- Does NOT increase EPL – still limited by mirror R
- Requires larger mirrors than axial configurations


What is cavity-enhanced spectroscopy?



Alternate approach: Ringdown spectroscopy



Alternate approach: Ringdown spectroscopy



$$\tau_0 = \frac{nL}{c(1 - R)}$$

$$\tau = \frac{nL}{c(1 - R + \alpha L)}$$

$$\alpha = c \left(\frac{1}{\tau} - \frac{1}{\tau_0} \right)$$

Comparison of CEAS techniques

NICE-OHMS

- Best laboratory performance, not field-ready

OA-ICOS

- Good laboratory performance, best field performance
- Large mirrors → large sample volume
- Significant performance loss if cavity is shortened

CRDS

- Best laboratory performance, good field performance
- Small mirrors → small sample volume
- Can be short (10 cm), but needs piezo actuators to account for large FSR

Design tradeoffs: Mirror reflectivity

- Higher reflectivity mirrors give longer effective path length, but make coupling more challenging
- Available reflectivity varies by wavelength – best at telecom bands (>99.999%), worst in UV / thermal IR (>99.97%)

$$\frac{12 \text{ cm}}{1 - 0.9997} = 400 \text{ m}$$

$$\frac{12 \text{ cm}}{1 - 0.99999} = 12,000 \text{ m}$$

Design tradeoffs: Cavity length

- Challenging to measure $\tau < 1 \mu\text{s}$; for 99.97% mirrors, minimum cavity length $\approx 10 \text{ cm}$
- Longer cavity more sensitive, but larger volume and may be difficult to fit in small spacecraft
- Longer cavity gives better FSR – could do without piezos in a large instrument, not practical for compact instruments

Design tradeoffs: Triggering method

Acousto-optic modulator:

- Fast response, does not perturb laser
- High power, water cooled


Semiconductor optical amplifier:

- Fast response, does not perturb laser
- Limited wavelength availability


Laser current modulation:

- Slower response, laser requires re-stabilization period
- No extra equipment, available for all semiconductor lasers

Prototype Water Analyzer Overview




Cell design




- 12 cm mirror separation
- 630 m EPL
- Mirrors masked to 5mm \varnothing
- 2.35 cc probed volume

Cell design: Plano-concave cavity



Cell design



FPGA system

FPGA

- Modulates laser
- Modulates piezo
- Interfaces with ADC
- Detects & triggers RDs
- Records RDs & metadata
- Passes data to CPU


CPU

- Sets laser scan parameters
- “Locks” piezo range
- Fits RDs to exponential function
- Controls fringe cancellation system
- Writes spectrum & metadata to file

Fringe cancellation


$$L = \frac{n\lambda}{2}$$

$$L = \frac{n\lambda}{2}$$




$$L = \frac{n\lambda}{2} - \frac{\lambda}{4}$$


$$L = \frac{n\lambda}{2} + \frac{\lambda}{2}$$




Fringe cancellation




Fringe cancellation




Spectrometer performance




Spectrometer performance



Spectrometer performance



Spectrometer performance



Spectrometer performance

Scan Parameter	Value
Cell Pressure	1 mbar
H ₂ O Partial Pressure	45 μ bar
Temperature	~294 K
Scan Time	50 s
Laser Power	10 mW

Isotope Ratio	Precision (‰ , 1σ)
$\delta^{18}\text{O}$	± 0.024
$\delta^{17}\text{O}$	± 0.11

What can we improve in the future?

- Demonstrate performance at other wavelengths
 - Currently working on methane measurements at $3.27\mu\text{m}$
- Demonstrate reliability in field settings – UAV measurements
 - Need to design and build second-generation cell, miniaturized electronics
- Find source for semiconductor optical amplifiers at non-telecom wavelengths

Acknowledgements

Mini-TLS Team:

- Chris Webster
- Jordana Blacksberg
- Didier Keymeulen
- Greg Flesch
- Erik Alerstam



Jet Propulsion Laboratory
California Institute of Technology