

# 67<sup>th</sup> Denver X-ray Conference

August 6<sup>th</sup> – 10<sup>th</sup>, 2018

Westminster, Colorado, USA

Evaluating perspectives from past missions to shape future investigations using the Mars 2020 Planetary Instrument for X-ray Lithochemistry

**C. M. Heirwegh<sup>1\*</sup>**

B. C. Clark<sup>2</sup>, A. C. Allwood<sup>1</sup>, D. T. Flannery<sup>1</sup>, J. A. Hurowitz<sup>3</sup>  
and W. T. Elam<sup>4</sup>

1) Jet Propulsion Laboratory, California Institute of Technology, 4800 Oak Grove Drive, Pasadena CA 91109

2) Space Science Institute, 4750 Walnut Street, Suite 2015, Boulder, CO 80301

3) Department of Geosciences, Stony Brook University, 100 Nicolls Road, Stony Brook, NY 11794

4) Applied Physics Laboratory, University of Washington, Seattle, WA 98105

\* - presenter



# XRF in Planetary missions



## NASA Contact science missions Mars employing XRF

- **Viking 2** (1975)
- **Pathfinder** (1996) – Sojourner
- **MER** (2003) – Mars Exploration Rover – Spirit and Opportunity rover
- **MSL** (2012) – Mars Science Laboratory – Curiosity rover
- **Mars 2020** – Planetary Instrument for X-ray Lithochemistry (PIXL)

$^{209}\text{Cd}$ ,  $^{55}\text{Fe}$

APXS

APXS

APXS

Rh X-ray tube

## NASA Contact science missions employing XRF

- **Venera 13/14, VEGA 2** (1983/84) – Venus
- **PHILAE Lander** (2014) – comet 67P/Churyumov-Gerasimenko
- **VIXL** – Venus Instrument for X-ray Lithochemistry Rh X-ray tube – Venus
- **Moondiver** (20??) – The Moon

$^{238}\text{Pu}$ ,  $^{55}\text{Fe}$

APXS


APXS?

APXS?


# Viking Lander program

PLANETARY INSTRUMENT FOR X-RAY LITHOCHEMISTRY

- Landers - Viking 1 (Chryse Planitia) and Viking 2 (Utopia Planitia)
- Launched 1975
- First Mars landers
- Spectrometer:  $^{55}\text{Fe}$  (5.9 keV) and  $^{109}\text{Cd}$  (22 keV) radioisotopes, proportional counter
- Rock chemistry, atmosphere, signs of life, images



**Figure:** Viking lander model. Source: [https://www.nasa.gov/sites/default/files/images/585497main\\_PIA09703\\_full.jpg](https://www.nasa.gov/sites/default/files/images/585497main_PIA09703_full.jpg)



**Figure:** Schematic of Viking lander sample analysis chamber. Source: [https://www.nasa.gov/sites/default/files/images/585497main\\_PIA09703\\_full.jpg](https://www.nasa.gov/sites/default/files/images/585497main_PIA09703_full.jpg)

# Viking 2 Lander - findings

PLANETARY INSTRUMENT FOR X-RAY LITHOCHEMISTRY



- Spectrometer resolution – 1000 eV @ 5.9 keV FWHM
- Weathered basaltic lava
- High **Fe** and Al, Si, K, Ca, Ti
- Very high **S** and **Cl** – ocean precipitates
- Rb, Sr, Y and Zr traces
- Soils same on both sides of planet (wind storm distribution)

## Deliverable

- First in-depth look at surface chemistry
- Co-validation of later data





Figure: Overlay of Viking measured soil spectrum with example APXS spectrum.




# Mars Pathfinder Mission

PLANETARY INSTRUMENT FOR X-RAY LITHOCHEMISTRY

- Rover “Sojourner” – launched 1996
- Landing Site: Ares Vallis
- First rover and airbag landing
- APXS – Alpha-proton X-ray Spectrometer
  - ~250 eV @ 6.4 keV



**Figure:** Pathfinder APXS spectrum of A-15 dark soil at Mermaid Dune site.  
Source: Rieder et al. 1997b Science **278** 1771-1774.




**Figure:** View of Sojourner rover source:  
<https://www.nasa.gov/specials/pathfinder20/>

## Discoveries

- Evidence of lakes and running water – Pebbles
- Mars Radius
- Atmospheric data
- Bulk Mars soil composition – similar to Viking. [S] and [Cl] differ.

**Figure:** Pathfinder APXS mounted on rover. Source: Rieder et al. 1997a J. Geophys. Res. E **102** 4027-4044.




# Mars Exploration Rover (MER)


PLANETARY INSTRUMENT FOR X-RAY LITHOCHEMISTRY

- Spirit and Opportunity Rovers – launched 2003
- Landing Sites: Columbia Memorial Station – Spirit, Meridiani Planum – Opportunity
- Opportunity – Active (2004 – present)
- Past water activity, rock chemistry, habitable regions
- APXS resolution – 160 eV @ 5.9 keV
- E cutoff – 16 keV

**Figure:** View from Cape Tribulation on rim of Endeavour Crater, Jan. 22, 2015. Source: <https://www.jpl.nasa.gov/spaceimages/details.php?id=pia19109>



**Figure:** MER rover in high-bay at JPL.  
Source:  
<https://mars.nasa.gov/mer/gallery/spaceraircraft/images/liverover.jpg>



## Discoveries & Comments

- Sulfates – hydrated [H<sub>2</sub>O] – Compton-Rayleigh
- Gusev
  - Outside rock coatings and cracks –water flow
  - Nickel – meteorites
- Columbia Hills
  - P, S, Cl, Br enrichment – aqueous alteration
- Clays, pure silica, carbonates – water exposure.
- “blueberries” – contains hematite

**Figure:** MER APXS spectrum of high-sulfate soil. Source: Campbell JL et al. (2008) J. Geophys. Res. 113 E06S11.

# Mars Science Laboratory (MSL)

PLANETARY INSTRUMENT FOR X-RAY LITHOCHEMISTRY



**Figure:** Artistic rendering of rover descent via sky crane. Source: [https://www.nasa.gov/mission\\_pages/msl/multimedia/gallery/pia14839.html](https://www.nasa.gov/mission_pages/msl/multimedia/gallery/pia14839.html)

Curiosity Rover – Launched 2011 – Gale Crater  
Mission: Habitability and humans to Mars

- Organic compounds
- Biosignatures
- Surface radiation
- Geological history

Figure: MSL Rover self image taken on Mount Sharp. Source: <https://photojournal.jpl.nasa.gov/catalog/PIA19808>.



Findings: Mars could have supported life

- Volcanic rocks - weathering
- Organic element building blocks
- Methane found
- Strong further evidence of water
- Cold, wet conditions prevailed
- Martian dust – high S & Cl




# Mars Science Laboratory (MSL) – cont'd


PLANETARY INSTRUMENT FOR X-RAY LITHOCHEMISTRY

- APXS – scientific methods development – MSL scientists
  - Relief effects
  - Oversampling - chemical lithology
  - Mineral phase effects (Na, Mg, Al, Si)
  - X-ray diffraction
  - C/R –ratios – “invisible” elements
  - Characterizing detector non-linearity – fits

- APXS properties:
  - Resolution (140 – 180 eV)
  - Full Pu L X-ray line energy range – improved C/R ratios



**Figure:** MSL APXS spectrum of a Martian soil. Source: Campbell JL et al. (2016) Nucl. Instrum. Meth. B 383 (2016) 143 – 151.



**Figure:** schematic of MSL APXS. Source: Campbell JL et al. (2013) Nucl. Instrum. Meth. B. **302**, 24 – 31..



# Perspectives

PLANETARY INSTRUMENT FOR X-RAY LITHOCHEMISTRY



- APXS method – trustworthy, heritage reliability
- APXS techniques – becoming advanced
- Geology and planetary chemistry – ongoing
- Habitability and biosignatures – more searching to be had
- Sub-mm scale probing – uncover individual grain chemistry

# Mars 2020 Mission Overview

PLANETARY INSTRUMENT FOR X-RAY LITHOCHEMISTRY


## Mars 2020 Science Objectives

- **Habitability** – past relevant environments
- **Biosignatures** – preservation potential
- **Sample caching** – potential future return
- **Prepare for humans** – *in situ* resources

Mission life: 1.5 Mars years  
(1005 earth days)

Launch: July, 2020


Land: February, 2021



## Instrument Complement:

- Mastcam-Z and Supercam for panoramic/stereo imaging and chemical analysis (LIBS)
- MEDA for weather
- RIMFAX ground penetrating radar
- MOXIE technology experiment to produce Oxygen from CO<sub>2</sub>
- SHERLOC and WATSON for UV Raman and high resolution imaging
- **PIXL micro-XRF spectrometer and Micro Context Camera (MCC) imager**

## Mars 2020



Total cameras: 23

Engineering cameras: 9

Science cameras: 7

Entry, descent and landing cameras: 7

# Planetary Instrument for X-ray Lithochemistry (PIXL) – Overview



PLANETARY INSTRUMENT FOR X-RAY LITHOCHEMISTRY

## Sensor Head





Figure (left): PIXL hexapod sensor head showing struts for (x,y,z) adjustments – source: L. Wade (JPL)

## Active Hexapod System (6 Active Struts)

## Key components:

- Micro-focus XRF
  - maps, lines, grids
- Optical fiducial subsystem (OFS)
  - exact (x,y,z) rock position
  - image co-registration
  - thermal drifting



| Parameter                | Performance Characteristics |
|--------------------------|-----------------------------|
| Excitation voltage       | 28 kV                       |
| Energy resolution        | 125 eV @ 5.9 keV            |
| Beam diameter (FWHM)     | 0.12 mm @ 8.0 keV           |
| Field of View            | 25 x 30 mm                  |
| Step size, any direction | 0.06 mm                     |

Table: PIXL operation characteristics – source: L. Wade (JPL)

Figure (center): PIXL sensor head front face. Two top cylinders – coarse and fine laser point grids, bottom center – micro-context camera, side circles – two Ketek SDDs, middle - x-ray optic.  
– source: L. Wade (JPL)



# Core Science Team & Locations

PLANETARY INSTRUMENT FOR X-RAY LITHO CHEMISTRY



Selen Turkey  
QUT



Chris Heirwegh  
JPL



ABIGAIL ALLWOOD  
PRINCIPAL INVESTIGATOR

JPL




JOEL HUROWITZ  
DEPUTY PRINCIPAL  
INVESTIGATOR

STONYBROOK



YANG LIU  
INVESTIGATION  
SCIENTIST

JPL



JOHN LEIF  
JORGENSEN  
Co-  
INVESTIGATOR

DTU



TIM ELAM  
Co-  
INVESTIGATOR

APL



DAVID FLANNERY  
COLLABORATOR

QUT



ROBERT HODYSS  
Co-  
INVESTIGATOR

JPL




BEN CLARK  
Co-  
INVESTIGATOR

SSI



MICHAEL TICE  
Co-  
INVESTIGATOR

U TEXAS AUSTIN



JOHN  
GROTZINGER  
Co-  
INVESTIGATOR

CALTECH



SCOTT  
MCLENNAN  
Co-  
INVESTIGATOR

STONYBROOK





ALLAN TREIMAN  
Co-  
INVESTIGATOR

LPI



## BHVO2 spectra – JPL breadboard prototype



# Measurement types: Line Scans


Short experiment, minimal positioning demands, quick science survey




1 x 24mm line: ~30 mins  
3 x 24mm lines: ~1hr

**Target:** Abraded / Natural  
**Suits:** 1 D heterogeneity / fast survey  
**Z:** Best achievable on natural surface, 25mm on abraded/flat  
**TOD:** Day (typical)  
**Goal:** Summed integration time on each component sufficient to determine the major/minor elements present and their relative variations.

## 1 Line



## 3 Line




## Ops Templates

**Step size:** 0.12 mm (=beam diameter) along lines  
**Line length:** 24 mm  
**No. of lines:** 1-3  
**OFS Imaging:** Intervals <4 mm, as shown, ensuring continuous coverage of dense SLI pattern for distance determination and X-ray spot location

○ Full-sized Image    ● Reduced Image    .... Spectral Line


# Measurement Types: Grid

Moderate duration, minimal positioning demands, more detailed science survey



**Target:** Abraded / Natural  
**Suits:** 2D heterogeneity moderately complex texture scale mm+  
**Z:** Best achievable on natural surface, 25mm on abraded/flat  
**TOD:** Day (typical)  
**Goal:** Summed integration on each component sufficient to determine major/minor elements and relative variations. Quantitative analysis and trace element detection possible on dominant phases

## 24 x 24 grid



## Ops Templates


Step size 1 mm (~10 x beam diameter) along lines  
Grid size: 24 x 24 mm  
OFS Imaging: Intervals <4 mm, as shown

○ Full-sized Image    ● Reduced Image    .... Spectral Line




# Measurement Types: Map

Greatest science value, longest duration, most demanding positioning requirements




**Target:** Abraded  
**Suits:** 2D heterogeneity Highly complex, very fine texture, SHERLOC overlap  
**Z:** 25mm  
**TOD:** Overnight  
**Goal:** Summed integration time on each component sufficient for quantitative analysis and trace element detection.


5 x 15 mm



7 x 7 mm



2 x 2 mm



○ Full-sized Image    ● Reduced Image    .... Spectral Line

**Ops Templates:**

**Step size:**


0.12mm (beam diameter)

**Dimensions:**

2 x 2mm,  
5 x 15mm  
7 x 7 mm


# Tactical Data Products: Summed, quantified spectra

PLANETARY INSTRUMENT FOR X-RAY LITHOCHEMISTRY



Spectra are summed, quantified: All spectra (bulk quantification) and selected rock components. Rock components selected by visible features, or **process shown on next slide**

## 1. Summed spectrum plot with peak labels



## 2. Table of calculated abundances

| Oxide                          | Wt % -<br>Comp 1 | Wt % -<br>Comp 2 | Wt % -<br>Comp 3 |
|--------------------------------|------------------|------------------|------------------|
| Na <sub>2</sub> O              | 7.11             | 2.08             | 1.5              |
| MgO                            | 3.61             | 0.82             | 5.65             |
| Al <sub>2</sub> O <sub>3</sub> | 16.09            | 5.66             | 5.29             |
| SiO <sub>2</sub>               | 50.7             | 73.65            | 30.93            |
| P <sub>2</sub> O <sub>5</sub>  | 0.49             | 1.25             | 0.88             |
| SO <sub>3</sub>                | 2.48             | 4.8              | 12.18            |
| Cl                             | 0.87             | 0.29             | 1.66             |
| K <sub>2</sub> O               | 2.23             | 0.96             | 0.09             |
| CaO                            | 6.07             | 3.05             | 25.26            |
| TiO <sub>2</sub>               | 0.49             | 1.57             | 0.52             |
| MnO                            | 0.14             | 0.06             | 0.87             |
| FeO                            | 9.47             | 5.49             | 14.61            |


# PIQUANT processing of PIXL data

PLANETARY INSTRUMENT FOR X-RAY LITHOCHEMISTRY

## Investigation tool (PIQUANT)

- FP
- Optic transmission
- Calculates spectra
- Quantifies (Z)
- **Work in progress**

Figure: PIQUANT fitted spectrum of NIST 610 with 60+ trace elements.



PIQUANT Version 2.01.126

☐ Energy Calibration    ☐ Plot Spectrum    ☐ Calculate Primary Spectrum    ☐ Calculate Full Spectrum  
☐ Compare measured to calculated    ☒ Fit one standard with plot  
☐ Calibrate    ☐ Quantify    ☐ Evaluate    ☐ Map

Configuration file:

Standards input file:

Calibration file:

Spectrum file:

Element fit controls Fe\_KLMN [IFX]:

Plot file (optional):

map file:

Log file (optional, always appends) ☒ Existing ☐ New

Command line options:

Fitted elements


|              |                    |               |              |
|--------------|--------------------|---------------|--------------|
| Na 10.1447 % | Na_K int 3610.6    | coeff 0.5662  | re_c 2.8%    |
| Mg 0.0465 %  | Mg_K int 137.2     | coeff 1.1901  | re_c 32.1%   |
| Al 1.0339 %  | Al_K int 10100.0   | coeff 1.2474  | re_c 1.1%    |
| Si 32.8733 % | Si_K int 629276.1  | coeff 0.9139  | re_c 0.1%    |
| P 0.0343 %   | P_K int -2130.2    | coeff -2.6921 | not included |
| Cl 0.0470 %  | Cl_K int 1924.6    | coeff 0.4644  | re_c 5.0%    |
| K 0.0461 %   | K_K int 8732.4     | coeff 0.6610  | re_c 4.4%    |
| Ca 8.3770 %  | Ca_K int 3894494.2 | coeff 0.9855  | re_c 0.1%    |
| Sc 0.0441 %  | Sc_K int 7193.8    | coeff 0.2833  | re_c 6.6%    |
| Ti 0.0437 %  | Ti_K int 42469.2   | coeff 1.2476  | re_c 1.2%    |
| V 0.0442 %   | V_K int 44200.3    | coeff 1.0171  | re_c 1.0%    |
| Cr 0.0415 %  | Cr_K int 53023.4   | coeff 1.0002  | re_c 0.8%    |
| Mn 0.0457 %  | Mn_K int 57314.2   | coeff 0.9047  | re_c 1.3%    |

- Elam, WT et al. (2002) Rad. Phys. Chem. **63**, 121-128.
- Ebel, H. (1999) X-ray Spectrom. **28**, 255-266.

# PIXL optic – Transmission profile derivation

PLANETARY INSTRUMENT FOR X-RAY LITHOCHEMISTRY

## Optic efficiency profile:



- Accuracy → FP approach
- Derive
  - XRF on H<sub>2</sub>O
  - PIQUANT calculated (no tube)




Figure (top): PIQUANT ECF's from fitted spectrum of NIST 610 with 60+ trace elements. Note divergence from unity (1).




Figure (bottom): Measured and PIQUANT-calculated water spectra.



# Optic profile results

PLANETARY INSTRUMENT FOR X-RAY LITHOCHEMISTRY

- New vs former profile differences
- Mostly corrected ECF trend
  - High Z ECF divergence cause due to fitting issue





Figure (right): New (smoothed) optic efficiency profile compared against former.

Source: **Heirwegh et al. (2018) Powder Diffraction Journal, 33, 162-165.**

Figure (left): Improved NIST 610 PIQUANT ECF agreement with unity. High Z divergence due to issues with background fit to spectrum.

# PIXL Quantification Capabilities

PLANETARY INSTRUMENT FOR X-RAY LITHOCHEMISTRY




## Quantification results:

- Informed PIQUANT tasks
- 8 powder standards: QLO1, JSM2, AN-G, SARM5, SARM6, BX-N, SY-4, GYP-B
- Avg error ~10% or better
- Next round – 16 powders

| Element            | Na | Mg | Al | Si | S | K  | Ca | Ti | Mn | Fe |
|--------------------|----|----|----|----|---|----|----|----|----|----|
| weighted error (%) | 9  | 2  | 4  | 2  | 9 | 17 | 4  | 10 | 16 | 7  |

Table: conc. Weighted error of eight powder certificate values vs PIQUANT concentration estimates.



# Selecting Rock sub-components


PLANETARY INSTRUMENT FOR X-RAY LITHOCHEMISTRY




**Dresser Formation, Western Australia ~3.5 Ga Microbialites,**  
cross cutting hydrothermal veins, exhalative hydrothermal  
deposits and multiple modes of alteration

**Bulk Summed Spectrum**

**Element Maps (net intensity) for  
Fe, Si, S, Ba, Ca peaks**



The element maps reveal variations in composition among rock subcomponents. To constrain the elemental composition, lithology and mineralogy of the subcomponents, we must group and analyze spectra from that subcomponent.




# Selecting Rock sub-components


PLANETARY INSTRUMENT FOR X-RAY LITHOCHEMISTRY




T-SNE (t-distributed stochastic neighbor embedding) scatterplot populations...



...back-projected on a context image



... and sums of the manually selected spectral populations





# PIXL: texture-specific chemistry

PLANETARY INSTRUMENT FOR X-RAY LITHOCHEMISTRY

## accurate correlation of chemistry with visible textures



Uses a focused X-ray beam to measure the chemistry of rocks, soils

Spatial resolution (beam  $\varnothing$ ) at focus = 0.12mm

X-Y scanning


Z motion to maintain focus

optical fiducial system (OFS) for context imaging and determining X-ray beam position on target



# Potential biosignature and context analysis



PLANETARY INSTRUMENT FOR X-RAY LITHOCHEMISTRY





# Greenland elemental maps

PLANETARY INSTRUMENT FOR X-RAY LITHOCHEMISTRY



# Concluding remarks

PLANETARY INSTRUMENT FOR X-RAY LITHOCHEMISTRY



- PIXL and Mars 2020 – continued discovery of Mars
  - Geological chemistry
  - Provenance – early development
  - Habitability
  - Biosignatures
  - Sample caching
  - Preparation for humans
- PIXL team efforts:
  - Ready flight instrument
  - Scientific preparedness
    - PIQUANT software
    - Elemental quantification calibration – FP approach
    - Map data visualization (eg. t-SNE)
    - Bulk quantification challenges (heterogeneities, rough surface, round edges)



# Acknowledgments





PIXL scientific data was acquired at the Jet Propulsion Laboratory, California Institute of Technology in Pasadena, California. Support for both data acquisition and the preparation and delivery of this talk was made possible through a contract with the National Aeronautics and Space Administration.

We extend our appreciation to the Denver X-ray Conference organization committee for their hard work in making this conference possible.

**Thank you for listening.**

# Tactical Data Products – Position vs. intensity plots


PLANETARY INSTRUMENT FOR X-RAY LITHOCHEMISTRY



- Mosaic and register MCC images with WATSON image
- Plot PIXL scan line(s) on WATSON image
- Process spectra
- Create position versus intensity plots (correlated to Watson image) for elements with adequate SNR


# Tactical Data Products – Element Intensity Maps

PLANETARY INSTRUMENT FOR X-RAY LITHOCHEMISTRY




Mosaic and register MCC images with WATSON image → plot PIXL scan lines on WATSON image → Process spectra → create net intensity plots

Net Intensity Plot (correlated to WATSON image) for elements with adequate SNR



PIXL scan area on WATSON image



Net Intensity Plot