



# Millimeter-wave Chirality Spectrometer (ChiralSpec)

Shanshan Yu, Theodore Reck, John Pearson, Michael Malaska and Robert Hodyss

Jet Propulsion Laboratory, California Institute of Technology, USA

**Brooks Pate**


University of Virginia, Charlottesville, VA 22903, USA




Copyright 2018. All rights reserved.

# Overview of ChiralSpec


- **Instrument description:** a simple mm-wave spectrometer operated in two modes:
  - Chirality detection mode to determine which enantiomer is in excess
  - Survey mode as a traditional mm-wave spectrometer to characterize chemical composition
- **Chirality detection theory formulated in 2012 by Hirota** and experimentally demonstrated by groups in Harvard and Virginia at 2-18 GHz with large instruments
- **Goal of this effort:** experimental demonstration of chirality detection at 75-205 GHz where there is a path for instrument miniaturization.
- **Addresses planetary science objectives** by looking for **biosignature patterns** : chirality of amino acids; double-bond position in fatty acids; amino acids distribution, etc.
- **Complimentary to mass spectrometer for distinguishing isomers of the same molecular weight**, e.g.,  $C_7H_5N$ , 103 Da, interesting to Titan, where is N inserted to? What is arrangement of atoms?




Benzonitrile



2-Ethynylpyridine



3-Ethynylpyridine




4-Ethynylpyridine

# ChiralSpec's core: rotational spectroscopy




- Non-invasive approach
- Extraordinary capability of distinguishing isomers via high resolving power




Multiple separate lines from each isomer ensure secured identification

- Capable of analyzing mixture of gases without separating them first



Survey spectrum of Orion KL from Herschel/HIFI (red numbers label  $CH_3OH$ )

# What is Chirality?





- Enantiomers of a chiral molecule: non-superimposable mirror images
  - Referred in terms of right- and left-handedness
  - Identical in many properties making them challenging to distinguish
  - Indistinguishable by mass spectroscopy
- Can be used as a chemical 'biomarker' for search for life
  - Organisms on Earth use predominately homochiral molecules
  - Life uses only left-handed amino acids
  - Life used only right-handed sugars

# Comparative Technologies

Technology	Methodology	Sensitivity	Comments
ChiralSpec	Rotational spectroscopy Targets in gas phase	Gas: $10^9$ molecules/cm <sup>3</sup> with 1L volume  Solid: ppm level with 0.2mg sample	Not require derivatization Insensitive to non-polar molecules
GC-MS	Gas chromatography Targets in gas phase	ppm	Require derivatization
MILA	Laser induced fluorescence Targets in aqueous phase	Liquid: 100 pM ( $10^{-10}$ mol/L)	


- Gas chromatography (GC), implemented in Rosetta/COSAC and Curiosity/SAM
- Microfluidic Life Analyzer (MILA)
  - Capillary electrophoresis-based laser induced fluorescence (CE-LIF)
  - Developed with PICASSO-13 (PI Peter Willis of JPL)
  - With focus on amino and carboxylic acids
  - Measurements involve liquid extraction of samples, fluorescent tagging, introduction of two chiral recognition agents

# Principle of chirality detection by ChiralSpec





- ✓ Exploiting **opposite signs of dipole moment products** (Image : Shubert et al 2016)
- ✓ Simultaneously using a cycle of three transitions involving the  $a$ ,  $b$  and  $c$  dipoles
- ✓ Generating a time-domain chiral emission signal via exciting species with two pulses
- ✓ Enantiomeric pair's chiral signals: same frequency but **180° phase difference**

# Schematic of ChiralSpec instrument



- All three waves capable of tuning frequencies in wide ranges
  - Capable of tuning to molecular fingerprint frequencies
  - Capable of **analyzing mixture of gases without separating them first**
- Chirality detection mode with all three waves on
  - Enantiomeric pair's chiral signals: same frequency but **180° phase difference**
  - Giving **difference abundance of enantiomeric pair**
- Survey mode with one wave on
  - Giving **sum abundance of enantiomeric pair**
  - Giving abundance for each of nonchiral molecules

# Benchmark Molecules

Molecule	Cycle	Pulse	Transition	Freq. (GHz)
Propylene Oxide <chem>(CH3CHCH2O)</chem> MW=58 g/mol $\mu_a, \mu_b, \mu_c = 0.95, 1.67, 0.56$ Debye 	1	Drive	$8_{4,4} \leftarrow 7_{3,4}$	182.856
		Twist	$8_{3,5} \leftarrow 8_{4,4}$	81.256
		Listen	$7_{3,4} \leftarrow 8_{3,5}$	101.599
	2	Drive	$15_{4,11} \leftarrow 14_{4,10}$	191.226
		Twist	$14_{5,10} \leftarrow 15_{4,11}$	88.069
		Listen	$14_{4,10} \leftarrow 14_{5,10}$	103.157
Alanine <chem>(CH3CHNH2COOH)</chem> MW=89 g/mol $\mu_a, \mu_b, \mu_c = 0.62, 1.33, 0.34$ Debye 	1	Drive	$21_{17,4} \leftarrow 20_{16,4}$	191.590
		Twist	$21_{16,5} \leftarrow 21_{17,4}$	76.855
		Listen	$20_{16,4} \leftarrow 21_{16,5}$	114.735
	2	Drive	$20_{19,1} \leftarrow 19_{18,1}$	195.558
		Twist	$20_{18,2} \leftarrow 20_{19,1}$	86.643
		Listen	$19_{18,1} \leftarrow 20_{18,2}$	108.915

- Propylene oxide is the first chiral molecule observed in ISM (McGuire et al 2016), is volatile and easier for sampling handling (vapor pressure of 450 Torr at 300K).
- We will use laser ablation to bring alanine into gas phase.





# Schedule and Milestones

Tasks	2017		2018				2019				2020			
	Q3	Q4	Q1	Q2	Q3	Q4	Q1	Q2	Q3	Q4	Q1	Q2	Q3	Q4
<b>T1 Design, fabricate, integrate &amp; test pulsing system</b>														
1.1 Drive pulse amplifier design,fabrication and test						◆ ~50mW achieved for the drive pulse power								
1.2 Phase-trigger circuit design, fabrication and test														
1.3 procurement														
1.4 Pulsing system integration and test										◆ 1 radian precision achieved for phase measurements				
<b>T2 Design, fabricate, integrate &amp; test resonator/coupler</b>														
										◆ Q ≥1000 achieved for resonators				
<b>T3 Advance mm-wave three-wave mixing technology</b>														
3.1 Measure chirality of propylene oxide without resonators													◆ Chirality measured for R&S sample	
<b>T4 Advance mm-wave cavity resonance technology</b>														
4.1 Measure chirality of propylene oxide with resonators						1000x chiral sensitivity increasement acheived								
<b>T5 Demonstrate ChiralSpec applicability to missions</b>														
5.1 Measure chirality of R-propylene oxide						Chiral signal observed for 10 <sup>-5</sup> mTorr R-sample							◆	
5.2 Measure chirality of mixture of R- and S-propylene oxide						Chiral signal observed for 10 <sup>-4</sup> mTorr sample with different R:S ratios							◆	
5.3 Set up laser ablation sample handling system														
5.4 Measure chirality of pure alanine ice														
5.5 Measure chirality of alanine-impregnated water ice						Chiral signal of alanine observed at 77K in alanine-impregnated water ice								◆



# Acknowledgements

---

The research was carried out at the Jet Propulsion Laboratory, California Institute of Technology, under Prime Contract NNN12AA01C issued through the PICASSO Program by the National Aeronautics and Space Administration.