



# EuroHEMS

Oct. 18, 2018

Cologne, Germany

## A Quadrupole Ion Trap for the Detection of Biomarkers at Icy Worlds

Josh Wiley<sup>1</sup>  
Anton Belousov<sup>1</sup>  
Sarah Waller<sup>1</sup>  
Dragan Nikolic<sup>1</sup>  
Richard Kidd<sup>1</sup>  
Stojan Madzunkov<sup>1</sup>  
Morgan Cable<sup>1</sup>  
Murray Darrach<sup>1</sup>

<sup>1</sup>*Jet Propulsion Laboratory, California  
Institute of Technology*


# Icy Worlds

- Ingredients for life (as we know it):
  - Water
  - Organic carbon
  - Energy source
  - Other elements: N, S, P
- 5 icy moons in our solar system that could potentially host life
- Need continued instrument development for lander and orbiter missions to search for signs of life

## Jupiter



Europa




Callisto




Ganymede

## Saturn



Titan




Enceladus

*Images: NASA*


# Enceladus

- Discovered by William Herschel in 1789 (named by son John in 1847)
- 500 km diameter
- One of the most likely icy worlds to host life
- Cryovolcanoes on the south pole eject ice, water vapor,  $H_2$ , NaCl, etc.
  - ~200 kg of material per second
  - ~500 km from the surface
- Widely accepted to have a large, salty, and subsurface ocean
- Over 100 geysers identified
- Ice particles make up Saturn's E-Ring





# Cassini: Enceladus' Ice Grain Analysis



Postberg, F.; *et al. Nature*, 558, p. 564, **2018**


Two mass analyzers were aboard Cassini

- Cosmic dust analyzer (CDA)
- Ion and Neutral Mass Spectrometer (INMS)

CDA:

- Time-of-Flight mass spectra of cations generated by high velocity (4-18 km/s) impacts of ice/dust grains on Rhodium target
- Mass range 1 – 200 Da (up to 8,000 Da)
- Mass resolution ( $m/\Delta m$ ): 20-50

## INMS Data



INMS:

- Detects ions and neutrals with quadrupole mass filter
- Mass range: 1 – 99 Da


Together, the CDA and INMS detected, water, salts (NaCl), benzene, methane, carbon dioxide and complex organics



# QIT-MS Development at JPL

2000


Quadrupole Ion Trap  
Mass Spectrometer  
(QIT MS)



Prior, Current or  
Funded Flight  
Instrument



Recent Work



2005

VCAM  
Flight  
GC




Gas Chromatography  
(GC)



VCAM  
Flight QIT MS


2010



MEMS GC  
Development




Vehicle Cabin  
Atmosphere Monitor  
(VCAM)




Wireless  
QIT MS

Noble Gas Isotopes  
0.1% accuracy




MARINE: ICCE  
Technology  
Maturation For  
MS




2015

HEOMD / AES  
Next Generation  
Life Support




Spacecraft  
Atmosphere  
Monitor




Cupid's Arrow  
Venus Probe



Deployment of  
Flyby MS on ISS




# QIT-MS on 4 New Frontiers Concepts




# Sampling Hypervelocity Ice with QIT-MS

## Ice Grain Sampling at Enceladus



- QIT-MS aboard spacecraft travelling several km/s
- Ice grains from Enceladus' ocean are sampled by QIT-MS

## Ice Grain Sampling in the Lab




- Liquid beam formed in vacuum
- Laser impinges on liquid beam and causes shock wave-induced anisotropic dispersion of small droplets
- Small droplets travel at a few km/s into QIT-MS

**Goal: Mimic ice grain size and relative velocity in the lab**




# Internal Ionization of Neutral Species Produced by Ice Grain Volatilization

## Step 1




## Step 2




## Step 3

Ions are mass analyzed by multi-stage mass spectrometry



# Modifications to QIT-MS for Low Vapor Pressure Molecule Analysis

Photo of Modified QIT-MS




1 cm<sup>2</sup> slot for ice grains

- Slot cut into QIT end cap to accept ice grains
- Trap coated with silicon-based SilcoTek coating to prevent analyte “sticking”
- Fatty/Amino solutions pumped directly into QIT-MS slot (QIT temperature: 125 °C)
- Adding a “reagent” allows chemical ionization (CI)
- Absence of a “reagent” results in electron ionization (EI)

Solution Injection in QIT-MS Diagram

Solution with  
Analyte (A)




EI – Low Pressure of R ( $<10^{-5}$  torr):

CI – High Pressure of R ( $>10^{-5}$  torr):


# Various Ionization Modes

## Electron Impact

- Solvent and analyte pumped into QIT
- After some time, solvent is pumped away but analyte remains
- Electron beam ionizes analyte




Solvent (S) +  
Analyte (A)




## Chemical Ionization

### Method 1:

- Solvent and analyte continuously pumped into QIT
- Solvent is used as reagent for CI




Solution:  
Reagent (R) +  
Analyte (A)




### Method 2:

- Solvent and analyte pumped into QIT
- After some time, solvent is pumped away but analyte remains
- Reagent gas for CI added in separate port





Solution:  
Solvent (S) +  
Analyte (A)





# EI vs. CI




## EI


- Causes significant fragmentation (“hard”)
  - Complicates mass spectra of mixtures
  - Excellent fingerprint for pure samples
- Parent peak: **Molecular cation**

## CI

- Can be “soft” or “hard” depending on reagents/samples
- CI of lauric acid with water is “hard”
  - Similar to EI spectrum
- Parent peak: **Protonated**



# “Softening” Chemical Ionization



- Difference in proton affinity determines extent of fragmentation


Water > Acetonitrile > Acetone

Molecule	PA (kJ/mol)	$\Delta$ PA (kJ/mol)
Lauric acid	815-825?	
Water	691	124 (1.3 eV)*
Acetonitrile	779.2	36 (0.4 eV)*
Acetone	812	3 (0.03 eV)*

\* Estimated values

- Mass spectra show less fragmentation as you work down from water to acetone
- Water CI shows more of the protonated parent ion than acetone
  - Likely due to reagent size
  - Should be able to further decrease fragmentation by using a reagent with a similar proton affinity as acetone but lower mass


# Pressure Dependence of CI



- CI is very sensitive to pressure
- Optimum CI with  $1 \times 10^{-5}$  torr acetone in QIT
- Sharp decrease in counts likely due to ion loss from collisions with acetone
  - A lighter CI reagent could possibly reduce signal loss




# Amino Acid CI in Water Mixtures




MW: 75 Da

PA: 886.5 kJ/mol



## CI: Water



## CI: Water/Acetone (99/1)




- Adding acetone to a majority water background gas “softens” ionization
- A water:acetone ratio of 99:1 provides best results
  - +/- one order of magnitude also OK
- Could select separate CI reagent for softer ionization of amino acids



Molecule	Proton Affinity (PA, kJ/mol)	$\Delta\text{PA w/Glycine}$ (kJ/mol)
Water	691	195.5 (2.0 eV)
Acetone	812	74.5 (0.77 eV)

# Analysis of Mixtures




Molecule	Mass (amu)
Glycine	75
Alanine	88
Serine	105
Methionine	149
Phenylalanine	165

- EI mass spectrum quickly becomes congested with fragment peaks (esp.  $m/z$  30 – 100)
- Can differentiate all 5 parent ions
  - Some fragmentation
- Can perform MS/MS for further structural confirmation


# MS/MS with Water Collision Gas

## CID agent: Water (with 1% Acetone)



- Minimal losses during isolation
- ~50-70% fragmentation efficiency
- Fragments at 133 and 104 amu
  - 133: Ammonia loss
  - 104: Formic acid loss
- Works well for all of the amino acids so far
  - Having some issues with fatty acids (lose water before activation)

# Design for Hypervelocity Ice Sampling



- Space venting to save space and power consumption
- Introducing water vapor with the acetone CI reagent vapor should help mitigate the effects of pressure fluctuations
- Potentially send two CI reagents
  - One for fatty & amino acids
  - One for amino acids

Molecule	PA (kJ/mol)
Lauric acid	815-825?
Glycine	886.5
Acetone	812
Ammonia	853.6
Chloromethylene	874.1

# Next Steps

- MS/MS to identify potential confounding isomers in a mixture




Glycine  
75 amu



Iso 1  
75 amu




Iso 2  
75 amu




Fragments

- Interface with hypervelocity ice gun




*Side-ionization-capable QITMS with hypervelocity input and gas delivery*



*Enclosure optimization for interface to hypervelocity ice gun.*






# Hypervelocity Ice Gun



# Generating Hypervelocity Ice Grains in the Lab

Absorption Spectrum of Water



IR Laser/Liquid Beam Interactions in the testbed generate neutral and charged particles with mass and velocity distributions analogous to the Enceladus Ice plume/QTIMS interactions.

# Conclusions

- Soft chemical ionization of amino/fatty acids was achieved in the QIT-MS in water-rich environments like those expected when sampling ice grains at hypervelocity
- A mixture of water:acetone (99:1) is capable of relatively soft ionization of both amino and fatty acids
  - Will continue to search for a lighter reagent with similar proton affinity as acetone
  - Could use a separate reagent for CI of amino acids with reduced fragmentation
  - Effective for mixtures
- The water/acetone mixture has been successfully used as a collision gas for MS/MS

# Acknowledgements

- Frank Postberg – Hypervelocity Ice Gun - Heidelberg
- Bernd Abel – Hypervelocity Ice Gun - Leipzig
- Other JPL team members:
  - All other authors, especially **Sarah Waller** and **Anton Belousov**
  - Jurij Simcic – Engineering/Design
  - Rembrandt Schaefer – RF Electronics
  - James Lambert – Hypervelocity ice (JPL)
  - Nick Tallarida – Hypervelocity ice (JPL)
  - Rob Hartsock – Hypervelocity ice (JPL)

