


Organic Capillary Electrophoresis Analysis System (OCEANS) Subsystem of the European Molecular Indicators of Life Investigation (EMILI)

Peter A. Willis¹, M. Fernanda Mora¹, Jessica S. Creamer¹, Aaron Noell¹, Florian Kehl¹, Konstantin Zamuruyev¹, Mauro Sergio Ferreira Santos¹, Tony Ricco², Richard Quinn², Jennifer Stern³, and William Brinckerhoff³

¹NASA Jet Propulsion Laboratory, Caltech, Pasadena, CA 91109;

²NASA Ames Research Center, Moffett Field, CA 94035;

³NASA Goddard Space Flight Center, Greenbelt MD 20771.




Outline

I. Architecting An Organic Chemical Analyzer for a Europa Lander Mission: EMILI

II. Key Developments at JPL


III. OCEANS



I. Architecting An Organic Chemical Analyzer for a Europa Lander Mission: EMILI

EMILI Philosophy

Deliver Maximum Europa Science by Casting the Widest Possible Nets in “Molecular Space”



Key & Driving Requirements:

- Must handle gases and liquids
- Must perform separations
- Must perform MS
- Must meet Europa Lander SDT measurement requirements

EMILI chooses **highest TRL hardware** with **lowest development risk** that is **capable of achieving science**


- **NASA Goddard** - MS(MOMA), GC(SAM), Mission Operations
- **NASA JPL** - Liquid Extraction, Electrophoresis
- **NASA Ames** - Spaceflight Microfluidics
- **Honeybee Robotics** - Sample Handling (SAM)

European Molecular Indicators of Life Investigation

EMILI = OCEANS + GAPS

OCEANS - Organic Capillary
Electrophoresis Analysis System

NASA JPL
NASA Ames





GAPS – Gas Analysis
Processing System

NASA Goddard

OCEANS uses Electrophoresis

- Simplest technique for performing liquid-based separation science on a mission
- Has low power, low mass, high sensitivity, high efficiency
- Charged molecules move at different speeds in the liquid
- **Various detection modes possible**



II. Key Developments at JPL

Most capable published method for amino acid analysis on life detection missions

by Jessica S. Creamer, Talk 107-3 Monday, June 24th, 2:00 – 2:15pm

**analytical
chemistry**


Article
pubs.acs.org/ac

Enhanced Resolution of Chiral Amino Acids with Capillary Electrophoresis for Biosignature Detection in Extraterrestrial Samples

Jessica S. Creamer,¹ Maria F. Mora, and Peter A. Willis*




LODs now down
to 1 nM



Peaks:
1. D-His
2. D-Leu
3. D-Val
4. L-His
5. L-Leu
6. D-Ser
7. GABA;
8. L-Val
9. D-Ala
10. L-Ser
11. β -Ala
12. L-Ala
13. Gly
*Dye side products

In Situ Chemical Analyzers in Development at JPL



LEGEND

Subcritical water extraction: Extracts organics from the sample

Sample processing unit: Characterizes sample and prepares it for analysis


Capillary/microchip electrophoresis: Separates sample into components

Laser-induced fluorescence: Detects labeled amino acids at nanomolar concentrations

Contactless conductivity: Detects any charged species, inorganic and organic ions

Electrospray ionization-Mass spectrometry: Detects a broad range of organic species

In Situ Chemical Analyzers in Development at JPL



LEGEND

Subcritical water extraction: Extracts organics from the sample

Sample processing unit: Characterizes sample and prepares it for analysis

Capillary/microchip electrophoresis: Separates sample into components

Laser-induced fluorescence: Detects labeled amino acids at nanomolar concentrations


Contactless conductivity: Detects any charged species, inorganic and organic ions

Electrospray ionization-Mass spectrometry: Detects a broad range of organic species

Chemical Laptop

Propulsion Laboratory
California Institute of Technology

First portable and
fully automated
ME-LIF instrument



. All rights reserved

Portable & Automated Extractor



Cap
Dispenser


Sample
Funnel

Extraction Cell
(Ball Aerospace
/ L. Beegle)

10 mm



Capping
Motor

Rotor




Biosignatures in the Atacama Desert

Do not need desalting
Do not need preconcentration
Just analyze fluid from extractor directly




Next: Validation in the Field



Jet Propulsion Laboratory
California Institute of Technology

Atacama Desert – September 2019




Extractor

Chemical
Laptop



III. OCEANS

In Situ Chemical Analyzers in Development at JPL



LEGEND

Subcritical water extraction: Extracts organics from the sample

Sample processing unit: Characterizes sample and prepares it for analysis

Capillary/microchip electrophoresis: Separates sample into components


Laser-induced fluorescence: Detects labeled amino acids at nanomolar concentrations

Contactless conductivity: Detects any charged species, inorganic and organic ions

Electrospray ionization-Mass spectrometry: Detects a broad range of organic species

In Situ Chemical Analyzers in Development at JPL

OCEANS



LEGEND

Subcritical water extraction: Extracts organics from the sample

Sample processing unit: Characterizes sample and prepares it for analysis


Capillary/microchip electrophoresis: Separates sample into components

Laser-induced fluorescence: Detects labeled amino acids at nanomolar concentrations

Contactless conductivity: Detects any charged species, inorganic and organic ions

Electrospray ionization-Mass spectrometry: Detects a broad range of organic species

Organic Capillary Electrophoresis Analysis System (OCEANS)



Coupling to mass spectrometry



Coupling to MS is essential in discovering “unknown unknowns” on life detection missions.

This leads us to select capillary approach.


A robust, spaceflight capable electrospray source that can be interfaced to a spaceflight compatible sample injector is the key missing piece of the puzzle.

CESI-MS

partnership with SCIEX (Brea, CA) Pioneer / World Leader in CE and CE-MS




- 1- Cyt
- 2- b-Ala
- 3- Ade
- 4- GABA
- 5- His
- 6- G-G
- 7- Gua
- 8- G-G-G
- 9- Cyt
- 10- Gly
- 11- Ado
- 12- Phe-Val
- 13- Ala
- 14- AIB
- 15- L-L-L
- 16- Val
- 17- iVal
- 18- Ser
- 19- Leu
- 20- Glu
- 21- Asp
- 22- Guo
- 23- Ura
- 24- Thd
- 25- Urd




OCEANS CE-C4D-ESI-MS


Electrophoresis system for multi-mode detection



fully
automated




precise, reproducible
sampling




OCEANS CE-C4D


System Validation




Conditions: 250 μM inorganics, 50 μM amino acids; BGE: 5 M acetic acid; injection volume: 8 nL; separation voltage: 20 kV; capillary length: 67 cm (effective: 47 cm); 50 μM I.D. capillary.



Reproducibility (10 runs)
RSD: <1% for migration time
RSD: 1-10% for peak height




OCEANS CE-C4D



Jet Propulsion Laboratory
California Institute of Technology


System Validation: Carry Over



- 1) NH_4^+
- 2) K^+
- 3) Ca^{2+}
- 4) Na^+
- 5) Mg^{2+}
- 6) Li^+
- 7) $\beta\text{-Ala}$
- 8) GABA
- 9) His
- 10) Gly
- 11) Ala
- 12) AIB
- 13) Val
- 14) Iva
- 15) Ser
- 16) Leu
- 17) Glu
- 18) Asp


Conditions: 250 μM inorganics, 50 μM amino acids; BGE: 5 M acetic acid; injection volume: 8 nL; separation voltage: 20 kV; capillary length: 67 cm (effective: 47 cm); 50 μM I.D. capillary.

OCEANS CE-MS



- portable MS (Waters QDA)
with SCIEX nanoelectrospray




- overlay of extracted ion chromatograms



Conditions: 50 μ M amino acids; BGE: 10% acetic acid; injection volume: 8 nL; separation voltage: 20 kV; pressure: 5 psi; capillary length: 90 cm; 50 μ M I.D.; spray voltage: 1.5 kV.

OCEANS Science Traceability

Goal	Objectives	Investigations
Search for evidence of life on Europa	A: Detect and characterize any organic molecular indicators of past or present life	A1: Determine the types, relative abundances, and enantiomeric ratios of amino acids (AAs)
		A2: Determine the abundances and patterns (i.e., population distributions) of organic compounds in the sampled material, with an emphasis on identifying potentially biogenic characteristics and molecular complexity.
	B: Characterize habitability	B1: Characterize the inorganic composition of the sample for sources of energy, chemically reactive species that could degrade biosignatures, or ions that could aid in elucidation of sample origin



Advancing Technology Readiness Level (TRL)


Key TRL Issues [1]:

- end-to-end operational testing,
- system integration,
- component selection and validation,
- sterilization,
- environmental non-operational testing,
particularly radiation [2]


[1] Willis et al., Analytical and Bioanalytical Chemistry 2015, 407 (23), 6939-6963.

[2] Creamer et al., Electrophoresis 2018 (22), p. 2864-2871. doi: 10.1002/elps.201800274.

Fieldwork



Jet Propulsion Laboratory
California Institute of Technology




Parting Message

EMILI maximizes the science return of a Europa Lander mission, *and our chances of identifying life, should it be present.*

OCEANS has been specifically designed for this purpose using NASA systems engineering processes.


The guiding principle behind our development strategy is the minimization of risk.

Acknowledgements



Jet Propulsion Laboratory
California Institute of Technology

- JPL's Chemical Analysis and Life Detection Group (pictured here)
- Funding from NASA/JPL programs: ICEE-2, JNEXT, PICASSO, NAI, PSTAR, NPP, SBIR, R&TD, TEFIM, Europa Lander SDT
- SCIEX
- The work done in this presentation was done at the Jet Propulsion Laboratory, California Institute of Technology, under a contract with the National Aeronautics and Space Administration.





Thank you

Recent Key Publications

“Stability of Reagents used for chiral amino acid analysis during spaceflight missions in high-radiation environments.” Creamer, J., Mora, M.F., and Willis, P.A. Electrophoresis 2018 (22), p. 2864-2871. doi: 10.1002/elps.201800274.

“Enhanced Resolution of Chiral Amino Acids with Capillary Electrophoresis for Biosignature Detection in Extraterrestrial Samples.” Creamer, J. S., Mora, M. F., & Willis, P. A., Analytical Chemistry 2017, 89(2), 1329–1337.

“Extraction of amino acids from aerogel for analysis by capillary electrophoresis. Implications for a mission concept to Enceladus’ Plume.” Mora, M. F., Jones, S. M., Creamer, J. and Willis, P. A. (2017), Electrophoresis 2017. doi:10.1002/elps.201700323

“Methods and Processes for Chiral Amino Acid Biosignature Analysis on Spaceflight Missions”, Jessica Creamer; Maria F. Mora; and Peter A. Willis, NASA NTR #50368, 2016.

“The Chemical Laptop: a Portable, Fully Automated Microchip Electrophoresis Instrument for Future in situ Planetary Investigations”, Peter A. Willis, Maria F. Mora, Florian Kehl, and Nathan Bramall, NASA NTR #50355, 2016.

“Implementation of microchip electrophoresis instrumentation for future spaceflight missions.” Willis. P. A, Creamer, J. S., Maria F. Mora. Analytical and Bioanalytical Chemistry 2015, 407 (23), 6939-6963.

Most Important Older Publications



Jet Propulsion Laboratory
California Institute of Technology

Cable et al., Microchip nonaqueous capillary electrophoresis of saturated fatty acids using a new fluorescent dye. *Analytical Methods* 2014, 6 (24), 9532-9535.

Cable et al., Identification of primary amines in Titan tholins using microchip nonaqueous capillary electrophoresis. *Earth Planet. Sci. Lett.* 2014, 403, 99-107.

Mora et al., Analysis of thiols by microchip capillary electrophoresis for in situ planetary investigations. *Electrophoresis* 2013, 34 (2), 309-316.

Cable et al., Low-Temperature Microchip Nonaqueous Capillary Electrophoresis of Aliphatic Primary Amines: Applications to Titan Chemistry. *Analytical Chemistry* 2013, 85 (2), 1124-1131.

Mora et al., Microchip capillary electrophoresis instrumentation for in situ analysis in the search for extraterrestrial life. *Electrophoresis* 2012, 33 (17), 2624-2638.

Mora et al., Toward Total Automation of Microfluidics for Extraterrestrial In Situ Analysis. *Analytical Chemistry* 2011, 83 (22), 8636-8641.

Willis et al., Monolithic photolithographically patterned Fluorocur (TM) PFPE membrane valves and pumps for in situ planetary exploration. *Lab on a Chip* 2008, 8 (7), 1024-1026.

Aubrey et al., The Urey instrument: An advanced in situ organic and oxidant detector for Mars exploration. *Astrobiology* 2008, 8 (3), 583-595.

Skelley et al., Organic amine biomarker detection in the Yungay region of the Atacama Desert with the Urey instrument. *Journal of Geophysical Research: Biogeosciences* 2007, 112 (G4).

© 2019. All rights reserved